

# Standard Operating Procedure

## CellTank 34

### 5 easy steps to single-use perfusion

1. Remove SUB from plastic bag
2. Connect sterile media bag and harvest bag to the SUB
3. Pump media into the SUB
4. Connect and calibrate the pre-installed classical SUS
5. Inoculate suspension or adherent mammalian cells and start cultivation

---

The following is recommendations, guide lines only. Originate 4<sup>th</sup> of January 2012, latest update December 2015. We suggest you develop your own Standard-Operation-Procedure (SOP).

---

#### Preparation of traditional **Re-Usable-Sensors (RUS)**:

- Install under sterile conditions on CellTank 34:
  - Required standard PG13.5x120 sensors; Glucose, Lactate, pCO<sub>2</sub> etc in appropriate *ports #10, 11, 12*
  - Biomass Hamilton/Fogale design only size PG13.5x120 sensor in *port #13* only and pre-amplifier with VP connector, M12 cabling to the instrumentation (not part of the CellTank kit) according to spec.
- Mount appropriate cabling to the installed RUS and calibrate appropriate the sensors.

#### Preparation with **Single-Use-Sensors (SUS)** in CellTank 34, which is: pH, dO<sub>2</sub>:

- Insure your have the appropriate cabling to the pre-installed classical pH SUS. Install and calibrate appropriate the sensor according to manufacturer spec.
- Install the VisiFerm PG13.5x120 sensor body (not part of the CellTank kit) without the; 1. Cap, 2.the tiny O-ring. 3. the large O-ring under the PG13.5 thread into the non-invasive well. Insure your have the appropriate cabling between the Hamilton SUS and you PCS. Install and calibrate the sensor according to the manufacturer spec.

#### Pre-cautions:

- DO NOT heat up the SUB until after media is added!
- Insure correct location of the SUB on a +25 watt Magnetic-Stirrer-Table (MST)
- DO NOT dry-run the integrated impeller with the MST!
- DO NOT stop the re-circulation pump at any time as this will stop the media flow through the matrix, scaffold containing the cells
- The product is designed for single-use only and will be damaged during treatment with pure ethanol or temperature above 50°C.

#### Media considerations

Serum containing media will create more foam compared to chemical defined (CD) media. Use of external foam collecting vessel size is depending on the serum content. 5 - 10% requires typically an external foam trap vessel for protection of the exhaust gas sterile filter. 2% serum containing media operates typically without external foam collecting vessel. Foam generation is heavily depending on the sparging gas volume and typically seen above 3 l/h.

#### Prep before start-up

1. Insert into VisiWell and calibrate the VisiFerm dO<sub>2</sub> sensor with air in the reservoir
2. Connect appropriate sized (>25 litre) media bag to port #2 and harvest bag to port #3

3. Add 1,000 ml media to the SUB under sterile conditions. The CellCore body volume is approx 280 ml.  $1,000 + 280 = \text{approx } 1,280 \text{ ml}$  level in vessel. Insure the sensors are fully covered with media. Remember the inoculation volume will add to the reservoir level!
4. Install the CellTank accurately at the centre on the clockwise rotating MST
5. Add 5-10 ml water or Glycerine to the non-invasive sensor port #8
6. Mount the thermocouple in the non-invasive port #8
7. Mount appropriate heating blanket around the vessel
8. Install the sky-support / filter support SS rod in the rod well on the cover
9. Expand the WC tube from port #16 vertically. Insure the hose clamp is locked. Insure the OD25 mm sterile filter is not blocked, mount the filter on sky-support.
10. Connect sparging gas supply to port #7 OD50 mm sterile filter, install optionally filter heater and secure the filter mechanically vertical
11. Connect head space flushing gas supply to port #6 OD50 mm sterile filter and secure the filter mechanically vertical
12. Initiate the MST and increase the rpm until you see the rotameter flotation device lifts inside the rotameter. Watch the inner CellCore is de-aired as gas and bubbles are exhausted through the tiny hole in the central rod.
13. Obtain stable condition at 0.8 litre/min flow (8 cm/min flux) on the re-circulation pump typically ranging 250-260 rpm.
14. Obtain stable conditions such as temperature
15. Calibrate the pre-installed pH Single-Use-Sensor according to the calibration chart found in the package.
16. Test three times for traction loss rpm in order to understand the MST potential. Increase with 5 rpm / second until the level is reached where the flotation device in the rotameter drops down and the SUB shakes. Make a note of the average traction loss rpm and keep 10% rpm distance below at any time.
17. Return to 0 rpm and increase to 0.8 litre/min flow typically 250-260 rpm, insure the media is flowing over the rotameter.
18. Open the hose clamp on the WC tube and insure the filter is not blocked.
19. Inoculate through port #1 or port #2 and start the cultivation in batch mode

CellTank's are pre-equipped with a 500 ml foam collection vessel, which may be installed hanging from the sky-support with filter heater on the exhaust sterile filter.

#### Gas control

1. Sparging gas ranging 3-6 l/h mixed gas 10-100% O<sub>2</sub> / 90-0% N<sub>2</sub>
2. Head space flushing gas range 4-8 l/h air, head space gas pressure >0.5 mBar for controlled pCO<sub>2</sub>, partial pressure in the media ~5 kPa.

#### Liquid level control

On Process-Control-Systems with no balancer but calibrated peristaltic pumps:

1. Use one pump to supplying fresh media (100% volume) to one of several *inlet port*. Harvest the product containing media (95% volume) from the bioreactor to suitable harvest storage through *harvest port*. Use one pump (10% volume) connected to *liquid level port* which then insures constant media level with return to harvest storage.
2. Use the pre-installed level sensor to control a peristaltic pump on/off or rpm operation. Remember that with SUS there is a need for an artificial ground to one of the electrodes taken from the PCS cabinet ground.

#### Inoculation of cells

High density – for the new user!!

Insure the cells are suspended carefully at high density in 50-100 ml media. Cell numbers to be inoculated depends heavily on the cell line and range for CHO cells  $0.1-2 \times 10^6 \text{ cells/ml/matrix}$ .

You may choose one of several routes for inoculation:

1. Into the reservoir through *port #5 or 6*
2. Into the volume between the impeller and the matrix through *port #2*

The impeller operating for a 0.8 l/min will distribute uniformly the cells inside the matrix.

Low density – for the experienced user!

A shaker flask with 1,000 ml media is used to grow up the cell line from frozen stage in an incubator at 37°C. Exchange the cap with a cap integrating a silicone hose and weld onto a suitable hose on the CellTank. Pump in the media in 10 minutes and set temperature on the PCS to 37°C. Start the re-circulation pump in the CellTank when +800 ml is transferred for 1 l/m flow.

### Separation of cells:

- Suspended cell line collection inside the matrix typically takes between 15 and 120 minutes. Hereafter the media is clear again and contains from 500 thousands down to thousands of cells/ml.
- Adherent cell lines separate by various collection phenomena inside the matrix and typically take between 5 and 50 minutes. Hereafter the media is again clear and contains from thousands down to only hundreds of cells/ml.

### Cell trapping control

The build-in impeller insures by rotation a constant re-circulation of media through the matrix.

- Suspension cell lines stay trapped inside the matrix cavities when media is constantly re-circulated at at least 0.8 l/m - avoid rapid media flow changes.
- Adherent cells stays trapped inside the matrix independent of the impeller is rotating or not – well, after seeding!

In any case media flow is needed in order to avoid gradients and volume determined by the cell mass, dO<sub>2</sub>, glucose and lactate content among others. Do not reduce the flux to a level below the requirements of the cells.




Figure – Cultivate in batch up to lactate / glucose decides perfusion start. Cultivate up to 50x10E+06 cells/ml/matrix and continue in perfusion. Reduce temperature with >5°C or more when lactate reaches 100% allowable and when glucose drops below acceptable levels.

### Batch cultivation

Start cultivation after inoculation in batch mode at >0.8 l/m re-circulation flow for 3-5 days and grow up to ~50x10E+06/ml cells/ml/matrix. Remember that correction factor between suspension cell to pF/cm is different from adherent cell lines.

The increasing high amount of cells will produce plenty of CO<sub>2</sub> and in order to reduce or reduce the addition of acid the head space must be flushed properly.

### **Suspension cell proliferation under perfusion mode operation**

Start the perfusion when glucose / lactate drops / increases beyond your acceptance. And at  $\sim 50 \times 10^6$  cells/ml/matrix reduce the SUB temperature from 37°C to a suitable level such as 32-30°C or lower in order to reduce proliferation.

Flux decrease dynamically according to increased cell density and range between 8 - 20 cm/min. Media replacement start from scratch and increase to as much as 30 matrix volumes per day depending on cell mass and cell type.

Expect a needed increase in impeller rpm each day in order to overcome the increased back pressure as to the increase cell density. Such as 10% increased rpm per day and an increase from start at 260 rpm to 500 rpm or more over 14 days operation. Be aware of the particular brand MST capability in use.

Maximum operational cell density is determined by the MST traction capability and varies between 500 and 750 rpm which are depending on brand! It is VERY important NOT to reach the traction loss rpm as the process will stop! Check traction guide lines.

DO NOT target  $1 \times 10^8$  cells/ml/matrix (unless you have such experience) before you reduce the SUB temperature.

### **Cell mass control**

There is several ways of controlling the maximum suitable bio mass in the matrix.

- Rapid increased impeller rpm (such as 50-100 rpm increase/second) will rip off/out high amount of suspension cells in the matrix to the reservoir. This high amount of cells in the reservoir is also known as “bleeding cells” from the SUB and may be taken out of the SUB through a suitable port. Cells in the reservoir will be trapped again after 15 and 120 minutes at suitable flux.
- Decrease the proliferation by altering the temperature to a lower level than 37°C on both adherent and suspension cell lines
- Decrease the proliferation by altering the glucose level on both adherent and suspension cell lines

### **Suspension cell perfusion cultivation**

Max cell density and the preferred method of keeping the cell density in control is the equilibrium for the production stage.



**Suggestions to the actual set-up**




Figure 4 – The illustrated set-up for CD media with four pumps controls the reservoir liquid level based on level sensed by conductivity via the level sensor.




Figure 5 – The illustrated set-up for serum containing media with five pumps determines the reservoir level as an alternative to a balancer or level control by conductivity. Port 9 leads foam to foam collection vessel which also then isolated the sterile exhaust filter from foam. Condensate foam is as a liquid re-directed to the vessel.

## Sampling

On three ports swabable valve with cover for sterile sampling are fitted.

#### **Data collection**

We suggest the following data is collected continuously. Further some suggested appropriate ranges for CHO cell lines:

- pH -
- dO<sub>2</sub> –
- biomass – 0-250 pF/cm (able to measure 700 pF/cm)
- Matrix sample, viable cell density
- Harvest sample, viable cell density
- Bottom sample, viable cell density
- Glucose – 10-20 mM
- Lactate – 5-50 mM
- Glutamine – 1.5-2-5 mM
- Glutamate – 0.5-1.5 mM
- Twice per day the re-circulation volume - 1-2.5 l/min
- Twice per day the WC – matrix pressure drop in mm Water Column – 100-500 mm
- Ammonia – 3-5 mM
- pCO<sub>2</sub> – partial CO<sub>2</sub> pressure 2-5 kPa in the reservoir controlled from head space
- Perfusion rate - reactor/volume – RV/day

PerfuseCell further recommend looking at the Application Notes and Posters  
END